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Terahertz RF Massager
ABS-WT01 / ABS-WT01A

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The technology values below are supplier-stated and must not be interpreted as medical claims. The user should follow intensity levels 1-10 and personal tolerance, not promotional temperature illustrations.

| | |
|------------------------|---|
| Model | ABS-WT01 / ABS-WT01A |
| Rated voltage | AC 220-240 V |
| Frequency | 50 Hz |
| Rated power | 400 W |
| Dimensions | 365 x 400 x 220 mm (L x W x H) |
| Net weight | Approx. 4 kg |
| Intensity levels | 1-10 |
| Timer | 5-30 minutes, adjustable in 5-minute steps |
| Ingress protection | IP20 - indoor dry use only |
| Electrical class | Class I |
| Test basis | EN 60335-1:2012/A16:2023; LVD Directive 2014/35/EU |
| Report | LCT-241102020332 |
| Terahertz band | Supplier-stated 3-3000 micrometer, with general frequency indication 0.1-10 THz |
| Electromagnetic module | Supplier-stated 1 MHz |
| Bioelectrical pulse | Supplier-stated 5-60 Hz |
| Pulse width | Supplier-stated approx. 200-400 microseconds |
| Magnetic field | Supplier-stated 2000 Gauss |

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- Heating belt: turns the heating belt on/off and selects temperature when connected.
- Red/Blue light: switches between red light, blue light and off.
- Time +/-: adjusts the operating time.
- Start/Pause: starts, pauses or ends the session.
- Intensity +/-: adjusts thermal/energetic intensity at levels 1-10.
- Plasma: turns the plasma function on/off.
- Time display: shows selected or remaining time.
- Intensity display: shows the selected level or "--" when foot contact is missing.
- Foot detection: alerts when one or both feet lose contact.
- Receiver: receives the signal from the remote control.

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- Check that the device is complete, clean and undamaged before use.
- Place the device on a dry, flat and stable surface with good ventilation.
- Do not use the device in bathrooms, wet rooms, near water or in high humidity.
- Remove metal objects such as ankle chains, jewelry and accessories that may come into contact with the device.
- The soles should be clean and slightly moist for good contact.
- Do not use the device if the skin on the feet is damaged, irritated, infected or has open wounds.

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- 1 Place the device level and stable. Make sure the power switch is off before connecting the cord.
- 2 Connect the power cord to the device and to a suitable outlet. Switch on the rear power switch.
- 3 In standby, the display usually shows 30 minutes and "--" in the intensity field.
- 4 Sit down and place both feet flat on the energy pedals without touching each other.
- 5 Press Start/Pause. When both feet have good contact, the device enters operating mode.
- 6 Select low intensity for first use and increase gradually only if comfortable.
- 7 Adjust the operating time as needed. Normal session length is 5-30 minutes.
- 8 When the time is finished, the functions switch off and the device returns to standby.
- 9 After use: switch off the device, unplug it, allow it to cool and clean if necessary.

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- Time and intensity: always start low and increase gradually. A long press on Start/Pause may switch off all functions.
- Red/blue light function: use with care and never look directly into the light source.
- Plasma function: the indicator flashes when active. Stop the session if you feel discomfort.
- Heating belt: optional accessory that must be connected only to the dedicated interface and used on dry skin. Do not wash or use the belt if it is damaged.

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- Pacemaker, cochlear implant or other electronic implant/active medical device
- Metal implants, stents or other implanted devices where use may be contraindicated
- Pregnancy, the first year after birth, infants, small children or persons unable to express discomfort
- Epilepsy, severe psychiatric illness, impaired consciousness or reduced ability to judge personal reactions
- High fever, infection, active skin infection, acute inflammation, bleeding tendency or severe unexplained symptoms
- Recent serious illness, heart failure, stroke, heart attack or other serious cardiovascular conditions
- Diabetes with reduced foot sensation, wounds, burns, swelling or pain in the affected area
- Light sensitivity, albinism or skin conditions where the red/blue light function should be avoided

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Nettoyage et stockage

- Switch off the device and unplug it after use.
- Allow the device to cool before cleaning.
- Wipe pedals and surfaces with a lightly damp, well-wrung cloth and then dry with a dry cloth.
- Do not use alcohol, solvents, abrasives or aggressive cleaning agents.
- Do not wash the heating belt and do not let it get wet.
- Store the device in a dry, clean place out of reach of children.

Dépannage

- Device will not switch on: check the outlet, power cord, rear power switch and fuse.
- Display shows "--": this normally means standby or missing foot contact. Lightly moisten the soles if needed.
- Device stops during use: the timer may have finished, foot contact may have been lost or a protection function may have activated.
- Remote control does not work: check battery, polarity, distance and clear line of sight to the receiver.
- Heating belt does not work: make sure the belt is connected correctly and remember it is controlled from the main panel.
- Unusual smell, smoke or strong heat: stop immediately, unplug and contact service.

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Conformité et documentation

- Test basis: EN 60335-1:2012/A16:2023
- Directive reference: LVD 2014/35/EU
- Report number: LCT-241102020332
- The product is not intended to diagnose, treat, cure or prevent disease.

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- Before use: read the manual, make sure the device is dry and undamaged, and sit down before starting.
- During use: keep both feet flat on separate pedals, start with low intensity and stop if uncomfortable.
- After use: switch off, unplug, let the device cool and wipe it with a lightly damp cloth.
- Do not use the device with pacemaker/implant, pregnancy, serious illness, skin damage, high fever or if you are unsure.